BAGNO STENDIFETTE ROTONDO

Codice: 40-300-002

CE IVD

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Rev. 01

PRECAUZIONI D'USO

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
Rev. 01

USE PRECAUTIONS

Before using the instrument, read carefully the instructions and warnings contained in this manual and keep it for further reference. They supply important indications regarding the functions and safety for installing, using and maintaining the instrument.

Bio-Optica Milano S.p.A. cannot be held responsible for any damage caused by improper or incorrect use and by the non-observance of any of the prescription provided in this manual and by the safety regulations in force.

1. After unpacking, make sure that the instrument is complete and not damaged by transport.
2. Before connecting the instrument to the power supply make sure that its rating corresponds to that of the power supply.
3. This instrument must only be used for the purpose for which it was designed, that is, as water bath for laboratory use. Any other use is to be considered improper and therefore hazardous.
4. The instrument must only be used by authorized and professionally qualified technician.
5. The electrical safety of this instrument can be guaranteed only if it is correctly connected to an efficient earth circuit as indicated by current electrical safety regulations. It is necessary to check this fundamental safety prerequisite, and if in doubt, ask to check the circuit. The instrument is provided with a power supply cable having 2 wires + ground tap that have to be connected to the power supply socket.
6. Do not remove the chassis or parts of it during operation. Switch off the instrument and disconnect the power supply cable before opening it. This operation must to be effected only by authorized and professionally qualified technician.
7. To eliminate instrument malfunctioning risks, do not work near strong magnetic fields and do not use transmitters such as cellular phones near the instrument. In case of serious malfunctioning switch off the instrument and contact the Technical Assistance Service.
8. All waste material, both infectious and radioactive, deriving from the appliance working cycle must be disposed in compliance with the regulation in force.

■ This appliance is marked from this symbol, in compliance with EU directive 2002/96/CE regarding electric and electronic appliances waste. This mean that the instrument, at the end of its useful life, must be collected separately from other refuse. The user must deliver it to the special differentiated refuse collection centres, that are predisposed by the public authority.

9. The contents of this manual is subject to change without further notice.
10. Please find enclosed the declaration of conformity.

11. Graphic symbols indicated on the label (positioned near the instrument's power supply socket):

Symbol for CATALOGUE NUMBER:

Symbol for SERIAL NUMBER:

Symbol for ALTERNATING CURRENT:

Symbol for FUSE:

Symbol for CONSULT THE INSTRUCTIONS:

Symbol for EC MARKING:

Symbol for IN VITRO DIAGNOSTIC-MEDICAL DEVICE:

Symbol for DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT:

Symbol for DATE OF MANUFACTURE:

Symbol for MANUFACTURER:

DESCRIPTION

Type of instrument: Thermostated bath for histology.

Dimensional features

Dimensions (∅ x H):	345x100 mm (external). 225x50 mm (basin).
Weight:	3,5 Kg.

Electrical connections

Power supply:	230V~ 50/60Hz.
Heating power:	200 Watt.

Other connections

Water connections:	Not necessary.
Fumes aspiration/filtration:	Not necessary.

Structural features

This multi-purpose bath is suitable for removing tissue in paraffin which come from microtome. The internal tank has been manufactured in black anodised aluminium making it easy to see the material contained inside. The external housing is made of oven painted aluminium coated with epoxy paint. The interspace has been insulated with thermal isolating to minimise heat loss and maintain uniform temperatures. Temperatures are regulated and controlled by a thermostat with fluid expansion sensor.

The instrument is provided with a safety thermostat integrated in the resistance that intervenes in case of malfunction or overtemperature above 100°C removing line voltage to the heater.

In the equipment it is included an anti-dust protection cover in anodized aluminium.

Technical features

Bath capacity:	1,5 lt.
Temperature:	Adjustable from +30°C to +80°C through analogic manual thermostat. Reading precision at +40°C, ±2°C.
Temperature taking system:	Through liquid expansion probe.
Heating power:	200W.
Protection rating:	IP54.

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Rev. 01

Controls

Switch ON/OFF	Switch for ignition and switching off of the instrument.
Warning light	It shows the working of the resistance.
Thermostat	It regulates the temperature of the sink.

Installation

Position the instrument on a level and stable working bench and connect it to the power supply socket (230V~ 50/60Hz) using the provided cable.

Important: Do not use any extension or adapter and do not modify the provided cable.

USE

Fill the basin at half height using a container. Switch on the instrument through the general switch ON/OFF.

Setting the temperature

Rotate the thermostat knob clockwise till the desired value is placed in correspondence of the notch. A green LED shows the ignition and switching off of the heating resistance.

Basin load/discharge

Switch off the appliance and unplug the cable from the socket before moving the instrument. Proceed with the replacement of the water in the sink.

Cleaning the metal casing and the basin

Before effecting the cleaning, switch off the instrument and unplug the cable from the socket.


Use only alcohol or non-aggressive detergents. Do not use abrasive products or acids in order to avoid ruining the varnished or plastic parts. The instrument does not need to be sterilized because is not expected the treatment of fresh samples but only of histological samples fixed and included in paraffin.

Maintenance

Not necessary.

Replacing the fuses

If the instrument doesn't start, check that the electric cable has been connected properly, there is current and the two fuses under the instrument's socket aren't burnt. If necessary, change them with fuses of identic value.

Fuses holder socket (N. 2 delayed fuses of 4 Ampere - T4AH250V).

Instructions for the replacement: Switch off the appliance and unplug the cable from the socket. Make a light pressure on the little carter covering the fuses (if necessary using a little screwdriver), change them, close the carter and check the instrument's ignition.


MAIN COMPONENTS

Code	Description	Quantity
37-WB001	sink	1
37-WB002	band	1
37-WB003	back	1
37-WB004	cover	1
37-WB006	seal for socket	1
37-POMWB100	knob	1
37-ETBIO.0	label	1
37-ET/HOT	"hot surface" label	1
37-PX6010063	plug	1
37-129-432	green waterproof switch	1
37-622	lamp	1
37-1023	rubber feets	4
37-1045	power supply cable	1
37-2090	thermostat	1
37-8010.0	resistance	1

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