PIASTRA CALDA PC400, PC800, PC810

Codici: 40-300-300, 40-300-301, 40-300-302

CE IVD

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


Code: 40-300-300, 40-300-301, 40-300-302

CE IVD

In the photo, from the right: Plate PC800, Plate PC810.

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Rev. 09

USE PRECAUTIONS

Before using the instrument, read carefully the instructions and warnings contained in this manual and keep it for further reference. They supply important indications regarding the functions and safety for installing, using and maintaining the instrument.

Bio-Optica Milano S.p.A. cannot be held responsible for any damage caused by improper or incorrect use and by the non-observance of any of the prescription provided in this manual and by the safety regulations in force.

1. After unpacking, make sure that the instrument is complete and not damaged by transport.
2. Before connecting the instrument to the power supply make sure that its rating corresponds to that of the power supply.
3. This instrument must only be used for the purpose for which it was designed, that is, as a hot plate for laboratory use. Any other use is to be considered improper and therefore hazardous.
4. The instrument must only be used by authorized and professionally qualified technician.
5. The electrical safety of this instrument can be guaranteed only if it is correctly connected to an efficient earth circuit as indicated by current electrical safety regulations. It is necessary to check this fundamental safety prerequisite, and if in doubt, ask to check the circuit. The instrument is provided with a power supply cable having 2 wires + ground tap that have to be connected to the power supply socket.
6. Do not remove the chassis or parts of it during operation. Switch off the instrument and disconnect the power supply cable before opening it. This operation must to be effected only by authorized and professionally qualified technician.
7. To eliminate instrument malfunctioning risks, do not work near strong magnetic fields and do not use transmitters such as cellular phones near the instrument. In case of serious malfunctioning switch off the instrument and contact the Technical Assistance Service.
8. All waste material, both infectious and radioactive, deriving from the appliance working cycle must be disposed in compliance with the regulation in force.

■ This appliance is marked from this symbol, in compliance with EU directive 2002/96/CE regarding electric and electronic appliances waste. This mean that the instrument, at the end of its useful life, must be collected separately from other refuse. The user must deliver it to the special differentiated refuse collection centres, that are predisposed by the public authority.

9. The contents of this manual is subject to change without further notice.
10. Please find enclosed the declaration of conformity.

11. Graphic symbols indicated on the label (positioned near the instrument's power supply socket):

Symbol for CATALOGUE NUMBER:

Symbol for SERIAL NUMBER:

Symbol for ALTERNATING CURRENT:

Symbol for FUSE:

Symbol for CONSULT THE INSTRUCTIONS:

Symbol for EC MARKING:

Symbol for IN VITRO DIAGNOSTIC-MEDICAL DEVICE:

Symbol for DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT:

Symbol for DATE OF MANUFACTURE:

Symbol for MANUFACTURER:

DESCRIPTION

Heated plate used for drying slides or for different laboratory uses.

Dimensional features

| | |
|---------------------|---|
| Dimensions (WxDxH): | 150x190x100 mm (PC400). 150x380x100 mm (PC800). 300x400x100 mm (PC810). |
| Weight: | 1,5 Kg (PC400). 2 Kg (PC800). 7 Kg (PC810). |

Electrical connections

| | |
|---------------|---|
| Power supply: | 230V~ 50/60Hz. |
| Power: | 85 Watt (PC400). 230 Watt (PC800). 400 Watt (PC810). |
| Fuses: | N. 2 fuses of 4 Ampere - 5x20 mm - T4AH250V. (PC800, PC810) N. 2 fuses of 1 Ampere - 5x20 mm - T1AH250V. (PC400) |

Structural features

- Thickness 15/10 painted sheet steel chassis.
 - Thickness 8 mm black anodized aluminium working plane.
-

Technical features

| | |
|---------------------------------|---|
| Working temperature adjustment: | From +20°C to +90°C through electronic thermostat with microprocessor. Reading precision $\pm 0,5^{\circ}\text{C}$. Control proportional of the temperature (PID). |
|---------------------------------|---|

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Rev. 09

Installation

Position the instrument on a level and stable working bench and connect it to the power supply socket (230V~ 50/60Hz) using the provided cable.

Important: Do not use any extension or adapter and do not modify the provided cable.

USE

Switch on the instrument by pressing the green main switch. The plate temperature will appear on display.


Setting the temperature

Press **BUTTON 1** for ca. one second in order to visualize the set temperature. By keeping **BUTTON 1** pressed, use **BUTTONS 2** or **3** to set the desired value. When **BUTTON 1** is released, the new value is stored.

BUTTON 4 is used by technical assistance service to set the parameters.

A red LED shows ignition and switching off of the heating resistance.

Frontal panel LTR-5 protection IP55

Changing the LAE thermostat parameters

Working parameters are preset according to reference standard values and can be modified.

In order to access the setup menu, press keys 1 and 4 for 5 seconds.

Select the parameter to modify using keys 2 and 3.

After choosing the interested parameter, press key 1 to visualize the value.

By keeping key 1 pressed, use keys 2 and 3 to set the desired value.

When key 1 is released, the new value is saved and the next parameter is visualized.

In order to exit from the setup, press key 4 or wait for 30 seconds.

| Parameter | Description | Value |
|-----------|-------------------------------------|---------------------|
| SPL | Minimum temperature that can be set | 20°C |
| SPH | Maximum temperature that can be set | 70°C |
| OS1 | Calibration of temperature probe* | between +01 and -02 |

*The standard deviation between the temperature shown by the probe and the real temperature is minimum ($\pm 1,5^\circ\text{C}$). In order to verify the real deviation, effect a measurement with a digital thermometer and compare the value with the temperature visualized on the thermoregulator.

If the temperature on the thermostat display is 40°C whereas the digital thermometer shows 38°C, subtract 2 from ADJ parameter in order to zero the deviation.

N.B. In order to guarantee the optimal working of the control proportional (PID) must not be modified the values relative to the following parameters: Pid, 1IT, 1DT, 1AR, 1CT.

Cleaning

Before effecting the cleaning, switch off the instrument and unplug the cable from the socket.


Chassis and plate: use only alcohol or non-aggressive detergents. Do not use abrasive products or acids in order to avoid ruining the varnished or plastic parts. The instrument does not need to be sterilized because is not expected the treatment of fresh samples but only of histological samples fixed and included in paraffin.

Maintenance

Effect the periodic cleaning.

Replacing the fuses

If the instrument doesn't start, check that the electric cable has been connected properly, there is current and the two fuses under the instrument's socket aren't burnt. If necessary, change them with fuses of identical value.


Fuses holder socket (N. 2 delayed fuses of 4 Ampere - T4AH250V).



Instructions for the replacement: Switch off the appliance and unplug the cable from the socket. Make a light pressure on the little carter covering the fuses (if necessary using a little screwdriver), change them, close the carter and check the instrument's ignition.

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