BAGNO STENDIFETTE WB1770

Codice: 40-300-000

CE IVD

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Rev. 11

PRECAUZIONI D'USO

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

Rev. 11

USE PRECAUTIONS

Before using the instrument, read carefully the instructions and warnings contained in this manual and keep it for further reference. They supply important indications regarding the functions and safety for installing, using and maintaining the instrument.

Bio-Optica Milano S.p.A. cannot be held responsible for any damage caused by improper or incorrect use and by the non-observance of any of the prescription provided in this manual and by the safety regulations in force.

1. After unpacking, make sure that the instrument is complete and not damaged by transport.
2. Before connecting the instrument to the power supply make sure that its rating corresponds to that of the power supply.
3. This instrument must only be used for the purpose for which it was designed, that is, as a water bath for laboratory use. Any other use is to be considered improper and therefore hazardous.
4. The instrument must only be used by authorized and professionally qualified technician.
5. The electrical safety of this instrument can be guaranteed only if it is correctly connected to an efficient earth circuit as indicated by current electrical safety regulations. It is necessary to check this fundamental safety prerequisite, and if in doubt, ask to check the circuit. The instrument is provided with a power supply cable having 2 wires + ground tap that have to be connected to the power supply socket.
6. Do not remove the chassis or parts of it during operation. Switch off the instrument and disconnect the power supply cable before opening it. This operation must to be effected only by authorized and professionally qualified technician.
7. To eliminate instrument malfunctioning risks, do not work near strong magnetic fields and do not use transmitters such as cellular phones near the instrument. In case of serious malfunctioning switch off the instrument and contact the Technical Assistance Service.
8. All waste material, both infectious and radioactive, deriving from the appliance working cycle must be disposed in compliance with the regulation in force.

■ This appliance is marked from this symbol, in compliance with EU directive 2002/96/CE regarding electric and electronic appliances waste. This mean that the instrument, at the end of its useful life, must be collected separately from other refuse. The user must deliver it to the special differentiated refuse collection centres, that are predisposed by the public authority.

9. The contents of this manual is subject to change without further notice.
10. Please find enclosed the declaration of conformity.

11. Graphic symbols indicated on the label (positioned near the instrument's power supply socket):

Symbol for CATALOGUE NUMBER:

Symbol for SERIAL NUMBER:

Symbol for ALTERNATING CURRENT:

Symbol for FUSE:

Symbol for CONSULT THE INSTRUCTIONS:

Symbol for EC MARKING:

Symbol for IN VITRO DIAGNOSTIC-MEDICAL DEVICE:

Symbol for DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT:

Symbol for DATE OF MANUFACTURE:

Symbol for MANUFACTURER:

DESCRIPTION

Thermostated bath used to collect and dry histological samples sections. Designed to be used together with the laboratory microtome, it maintains the water temperature at the desired level so that the sections can be handled properly.

Dimensional features

Dimensions (WxDxH):	350x365x155 mm.
Slide-dryer plate dimensions (WxD):	350x100 mm.
Weight:	8 Kg.

Electrical connections

Power supply:	230V~ 50/60Hz.
Power:	350 Watt.
Fuses:	N. 2 fuses of 3,15 Ampere - 5x20 mm - T3.15AH250V.

Structural features

- Thickness 15/10 painted sheet steel chassis.
 - NTC10K probe provided with articulated arm and safety device that cut off the power supply and the basin heating when the basin is removed.
 - Removable pirex basin (capacity ca. 1,2 litres) for a simple operation of filling and emptying.
 - A 6 Watt neon lamp illuminates the basin in order to facilitate the sections collecting.
 - Slide-dryer plate that can house up to 24 slides, situated on the back side.
-

Technical features

Working temperature adjustment (Water bath):	From +20°C to +70°C through electronic thermostat with microprocessor. Reading precision $\pm 1^\circ\text{C}$. Control proportional of the temperature (PID).
Temperature taking system:	Through NTC10K probe immersed directly in water.
Working parameters visualization and modification:	Through control panel with decimal digital display.
Working temperature adjustment (Slide-dryer plate):	From +20°C to +65°C through analogic manual thermostat. The full scale up to 90°C is not usable.

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Rev. 11

Installation


Position the instrument on a level and stable working bench and connect it to the power supply socket (230V~ 50/60Hz) using the provided cable.

Important: Do not use any extension or adapter and do not modify the provided cable.

USE

Instructions for the basin filling:

Pull the probe and rotate it clockwise (as shown in the picture).
Fill the basin with distilled water up to half level.
Immerse the probe in the water.

Switch on the instrument by pressing the green main switch. The water temperature will appear on the display.


Setting the Water Bath's temperature

Press **BUTTON 1** for ca. one second in order to visualize the set temperature. By keeping **BUTTON 1** pressed, use **BUTTONS 2** or **3** to set the desired value. We suggest to set it at ca. 40°C. When **BUTTON 1** is released, the new value is stored.

BUTTON 4 is used by technical assistance service to set the parameters.

A red LED shows ignition and switching off of the heating resistance.

Frontal panel LTR-5 protection IP55

Changing the LAE thermostat parameters

Working parameters are preset according to reference standard values and can be modified.

In order to access the setup menu, press keys 1 and 4 for 5 seconds.

Select the parameter to modify using keys 2 and 3.

After choosing the interested parameter, press key 1 to visualize the value.

By keeping key 1 pressed, use keys 2 and 3 to set the desired value.

When key 1 is released, the new value is saved and the next parameter is visualized.

In order to exit from the setup, press key 4 or wait for 30 seconds.

Parameter	Description	Value
SPL	Minimum temperature that can be set	20°C
SPH	Maximum temperature that can be set	70°C
OS1	Calibration of temperature probe*	between +01 and -02

*The standard deviation between the temperature shown by the probe and the real temperature is minimum ($\pm 1,5^\circ\text{C}$). In order to verify the real deviation, let the water contained in the basin reach the temperature of 40°C and measure the real temperature with a digital thermometer.

If the temperature on the thermostat display is 40°C whereas the digital thermometer shows 38°C, subtract 2 from ADJ parameter in order to zero the deviation.

N.B. In order to guarantee the optimal working of the control proportional (PID) must not be modified the values relative to the following parameters: Pid, 1IT, 1DT, 1AR, 1CT.

Setting the Slide-Dryer Plate's temperature

The plate heating starts by pressing the green main switch; to exclude it, rotate the thermostat knob anticlockwise up to end-run.

The instrument is tested and preset at a temperature of ca. 40°C; to modify it proceed as follows:

- 1) Switch on the instrument and wait ca. 15 minutes.
- 2) Measure the plate temperature with a thermometer or letting some slide dry.
- 3) Adjust the temperature rotating the black knob, making little movements with pauses of ca. 10 minutes.

Cleaning

Before effecting the cleaning, switch off the instrument and unplug the cable from the socket.

Pirex basin: use decalcifying products or acids for limescale, solvents for residual paraffin.


Chassis: use only alcohol or non-aggressive detergents. Do not use abrasive products or acids in order to avoid ruining the varnished or plastic parts. The instrument does not need to be sterilized because is not expected the treatment of fresh samples but only of histological samples fixed and included in paraffin.

Maintenance

Effect the periodic cleaning.


Opening the instrument

Switch off the instrument and unplug the cable from the socket. Remove the screws of the lower panel (the most recent models are assembled with special screws, type TORX with central hole). Replace the instrument on the feet and rotate the upper side towards left (as shown in the image).



Replacing the fuses

If the instrument doesn't start, check that the electric cable has been connected properly, there is current and the two fuses under the instrument's socket aren't burnt. If necessary, change them with fuses of identical value.

Fuses holder socket (N. 2 delayed fuses of 3,15 Ampere - T3.15AH250V).

Instructions for the replacement: Switch off the appliance and unplug the cable from the socket. Make a light pressure on the little carter covering the fuses (if necessary using a little screwdriver), change them, close the carter and check the instrument's ignition.

Replacing the lamp

Switch off the instrument and unplug the cable from the socket.

Remove the panel at the bottom of the instrument unscrewing the 8 fixing screws and lift the chassis rotating it on the left. Rotate the lamp of 90° and lift it. Change it with an equal model and rotate it of 90°. Close the instrument again.

Replacing LAE thermostat

Switch off the instrument and unplug the power supply cable from the socket, open the instrument. The thermostat is fixed with clips or blocking flasks. The connections are on the label positioned on the upper side of the thermostat.

TROUBLE SHOOTING

The instrument doesn't switch on

Check that the cable is fixed correctly, there is current in the electric system and the power supply is correct. Check the fuses working and change them, if necessary.

The instrument lights up but the keys doesn't respond to controls

It is possible that, after some time, the membrane covering the keys warps and blocks the keys. Replace the thermoregulator LAE.

The instrument lights up but the thermostat doesn't light up and the resistance doesn't heat

The "L"-shaped stainless steel arm containing the temperature probe in the working position close a microswitch that feeds the thermostat. Verify the microswitch working, the correct positioning and fixing of the internal disk.

An error signal (for example DL1) appears on the display

Verify the correct connection of the temperature probe.

The basin doesn't heat

The thermostat indicates with a red point the circuit closing and so the resistance power supply. If the red led doesn't light up, verify the thermostat parameters. If, on the contrary, it lights up check the power supply and the working of the heating resistance.

The slide-dryer upper plate doesn't heat

Verify the correct working of the analogic thermostat and the connection of the wiring. Check the power supply and the working of the heating resistance.

The lamp doesn't work

Replace the lamp or the starter.


MAIN COMPONENTS


CODE	DESCRIPTION	QUANTITY
40-300-050.0	pirex basin	1
37-AR09E2F6A	6 ampere wire filter	1
40-300-051.0	opaline white cover	1
37-MOL1770	steel spring	1
37-OPA220X30	white opaline	1
37-129-432	green waterproof switch	1
37-1004	50 watt resistance	1
37-1005	90° tube for feeler	1
37-1006/A	knob with pin	1
37-1006/B	threaded knob	1
37-1006/C	∅ 30x10 disk	1
37-1007	microswitch	1
37-1008	metal ring	1
37-1011	lamp holder	2
37-1013.0	reactor	1
37-1014	neon lamp	1
37-1023	rubber feets	4
37-1028	stirrup	1
37-1040/N	digital thermoregulator	1
37-SO1040/N	probe for thermoregulator	1
37-RELE1040-25	relay	1
37-1045	power supply cable	1
37-2090	thermostat	1
37-4090	polycarbonate panel	1
37-6010	mica resistance	1
37-1020/1	upper chassis	1
37-1020/2	black slide holder plate	1
37-1020/3	internal slide holder plate	1
37-1020/4	black pirex holder dish	1
37-1020/5	internal pirex holder dish	1
37-1020/6	lower panel	1

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