

Osteodec



In-vitro diagnostic medical device

EMDN : W01030799

BASIC UDI: 080339762W01030799Y5

IVD in **Classe A**, Reg. UE 2017/746

Catalog number	Unit size	UDI-DI
05-M03005	500 ml	08033976234416
05-03005Q	4 x 2,5 l	08033976232115

Packaging

05-03005Q

Primary container: white bottle in polyethylene terephthalate (PET). Useful capacity 2.5 liters. HDPE cap.

Tamper evident cap.

The polyethyleneterephthalate is a thermoplastic polymer of the polyester family. PET is an optimal oxygen, carbon dioxide and other gasses barrier. This material has a high resistance to ultraviolet radiation and an inertia toward the mainly chemical agents (solvents: xylene, limonene, liquid paraffines, alcohols, acids, bases etc.). It is biologically inert. It constitutes a good water and humidity barrier. It shows a great hardness and mechanical resistance.

The bottle has an optimal grip. The absence of the handles reduces space for storage. The anti-dropping cap permits a precise and clean use.

Secondary container: carton box.

05-M03005

Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 500 ml. HDPE cap. Tamper evident cap.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim

Product for the preparation of cyto-histological samples for optical microscopy.

Application

Decalcifier for needle biopsy, small bone fragments, calcified arteries, fleshy tissue pathologically mineralized.

Principle

This reagent is characterized by bisodic EDTA in acid buffer. Decalcification takes place through chelation in acid environment. Thanks to the co-existence of two modes of action, decalcification is quicker than usual EDTA procedure, preserving at the same time tissue structure.

- Fixation technique**
- 1) Volume ratio specimen/decalcifier 1 : 100
 - 2) Procedure time: 1 – 4 hours for needle biopsy and small fragments; 8 – 16 hours for 3-5 mm thick specimens.
 - 3) Post decalcifying procedure: reconditioning of the sample. Rinse with slowly running water for 1 hour or apply 3 changes of PBS pH 7.4, 20 minutes each.

Components

Components	CAS	CE	Index
E.D.T.A. disodium salt	60-00-4	2004494	-
Acid buffer			

Warning and precaution

The product is intended for professional laboratory use for healthcare professionals. Carefully read the information on the label (danger symbols, risk and safety phrases) and always consult the safety data sheet. Do not use if the primary container is damaged. In the event of a serious accident, we recommended that you immediately inform Bio-Optica Milano S.p.A and the competent authorities.

Storage

Store the preparation at room temperature. Keep the containers tightly closed.

Stability

After the first opening, the product is reusable until the expiry date, if correctly stored. Validity: 2 years.

Disposal

Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022