



Papanicolaou EA65

Cytoplasmic staining



In-vitro diagnostic medical device

EMDN : W01030708

BASIC UDI: 080339762W01030708X8

IVD in **Classe A**, Reg. UE 2017/746

Catalog number	Unit size	UDI-DI
05-12017	500 ml	08033976232771
05-12017/L	1 l	08033976232788

Packaging 05-12017
Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 500 ml. HDPE cap. Tamper evident cap.

05-12017/L
Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 1 l. HDPE cap. Tamper evident cap.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of: gynecological specimens, urine cytology, fine needle specimens, sputum and bronchial washings, to be examined by optical microscopy.

Application Cytoplasmic staining solution for Papanicolaou method.
For the execution of the staining method is required the use of reagents Papanicolaou Harris hematoxylin and Papanicolaou OG6.

Principle A highly selective blue nuclear stain, Harris' hematoxylin, is combined with EA65 polychromic mixture, a subtle cytoplasmic stain which differentiates cyanophil cells from eosinophil ones. The last ingredient is OG6 solution, which stains keratinized elements.

Method

- 1) Ethanol 95°, 2 minutes
- 2) Distilled water, 2 minutes
- 3) Harris Hematoxylin, 1 minute
- 4) Tap water, 5 minutes
- 5) Ethanol 95°, 15 seconds
- 6) OG 6, 2 minutes
- 7) Ethanol 95°, 15 seconds (twice)
- 8) EA 65, 5 minutes

- 9) Ethanol 95°, 15 seconds
- 10) Absolute Ethanol, 30 seconds (twice)
- 11) Xilene or Bio Clear, 2 minutes (twice)

Results

Nuclei : Blue-purple
 Cyanophil cytoplasm : Blue-green
 Eosinophil cytoplasm : Pink
 Keratinized cytoplasm : From pink to orange

Components

Components	CAS	CE	Index
Eosin Y Certified - CI 45380	17372-87-1	241-409-6	-
Light green - CI 42095	5141-20-8	225-906-5	-
Phosphotungstic acid	12501-23-4	-	-
Ethanol 95°	64-17-5	200-578-5	603-002-00-5

Warning and precaution

The product is intended for professional laboratory use for healthcare professionals. Carefully read the information on the label (danger symbols, risk and safety phrases) and always consult the safety data sheet. Do not use if the primary container is damaged. In the event of a serious accident, we recommended that you immediately inform Bio-Optica Milano S.p.A and the competent authorities.

Storage

Store the preparation at 15-30°C. Keep the containers tightly closed.

Stability

After the first opening, the product is reusable until the expiry date, if correctly stored. Validity: 2 years.

Disposal

Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022