

Killik - O.C.T. Compound
Embedding medium for cryostat



IVD In-vitro diagnostic medical device
EMDN : W01030799
BASIC UDI: 080339762W01030799Y5
IVD in **Classe A**, Reg. UE 2017/746

Catalog number	Unit size	UDI-DI
05-9801 (neutral)	4 x 100 ml	08033976233020

Packaging Primary container: neutral bottle in PE with resealable spout designed to an easy dispensing. Useful capacity 100 ml.
Secondary container: carton box.

Wear, water, alcohol and solvents resistant label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of cyto-histological samples for optical microscopy.

Application Inclusion matrix for the preparation of the tissue sample to the cut to the cryostat.

Method
1) Cover the sample with the embedding medium.
2) Freeze the sample.
3) After freezing, proceed with cutting.

Components	Components	CAS	CE	Index
	Water soluble polymers mixture	-	-	-

Functional properties

- The polymer's freezing temperature is optimal for sectioning at the cryostat.
- Water soluble
- Leaves no residues
- Adequate viscosity to support specimen
- Improves tissue adhesion to object holder
- Long-spout dispenser bottle, ideal for easy dispensing.

Warning and precaution
The product is intended for professional laboratory use for healthcare professionals. Carefully read the information on the label (danger symbols, risk and safety phrases) and always consult the safety data sheet. Do not use if the primary container is damaged. In the event of a serious accident, we recommended that you immediately inform Bio-Optica Milano S.p.A and the competent authorities.

Storage Store the preparation at room temperature. Keep the containers tightly closed.

Stability After the first opening, the product is reusable until the expiry date, if correctly stored.
Validity: 2 years.

Disposal Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022