

Distilled water – B-PRO 450



For Bio-Optica B-PRO 450 processor

IVD In-vitro diagnostic medical device
EMDN : W01030799
BASIC UDI: 080341202W01030799BM
IVD Class A, Reg. UE 2017/746

Catalog number	Unit size	UDI-DI
450003	5 l	08034120278492

Packaging Primary container: PE tank in neutral colour, capacity 5 liters, UN approved. Screwcap in PE with a seal, blue colour. Watertight.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Specifications		
	Liquid	clear, colorless
	Phosphates	≤ 0,1 ppm
	Ag	≤ 0,01 ppm
	B	≤ 0,01 ppm
	Bi	≤ 0,01 ppm
	Co	≤ 0,01 ppm
	Fe	≤ 0,01 ppm
	Li	≤ 0,01 ppm
	Mo	≤ 0,01 ppm
	Pb	≤ 0,01 ppm
	Ti	≤ 0,01 ppm
	Zr	≤ 0,01 ppm
	Nitrates	≤ 0,1 ppm
	As	≤ 0,01 ppm
	Ba	≤ 0,01 ppm
	Ca	≤ 0,01 ppm
	Cr	≤ 0,01 ppm
	Mg	≤ 0,01 ppm
	Na	≤ 0,01 ppm
	Sn	≤ 0,01 ppm
	V	≤ 0,01 ppm
	Si	≤ 0,01 ppm
	Chlorides	≤ 0,1 ppm
	Sulphates	≤ 0,1 ppm
	Au	≤ 0,01 ppm
	Cd	≤ 0,01 ppm
	Cu	≤ 0,01 ppm

K	≤ 0,01 ppm
Mn	≤ 0,01 ppm
Ni	≤ 0,01 ppm
Sr	≤ 0,01 ppm
Zn	≤ 0,01 ppm
pH	20°C 5.5 -7.5
Evaporation residue	≤ 1 ppm
Calcination residue	≤ 1 ppm
Conductivity at the time of production	≤ 0,1 μS/cm

Storage	Store in a cool place away from light. Keep the containers tightly closed.
Stability	After the first opening, the product is usable until the expiry date, if correctly stored. Product validity: 1 year.
Disposal	Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

REVISION N°	REASON	REVISION DATE
001	First issue	28/12/2022