

# **SAFETY BOX WITH FOAM INSERT FOR KLESSIDRA**



CODE	DESCRIPTION
05-900700	Safety box complete with foam with 21 places – 1 pcs
05-900800	Accessory: foam with 21 places – 1 pcs
05-900850	Accessory: foam with 5 places – 1 pcs
44-6754	Accessory: Safety container (460 x 240 x 230 mm) for code 05-900700

Distributor: Bio-Optica Milano S.p.A.

Date of emission: 30/08/2024

Rev. 01

### DESCRIPTION

Safety box for Klessidra complete with foam support with 21 places.

Made of polycarbonate. The cover has a hermetic silicon seal, to ensure a high degree of safety, even in case of accidental leak of liquids contained within, and it is provided with 4 safety locks preventing opening if the box drops down.

Complete with stainless steel handle. Biohazard logo printed on the lid.

### TECHNICAL FEATURES

Safety box	Dimensions	Height	180 mm
		Length	330 mm
		Depth	175 mm
	Material	Polycarbonate, colourless	
Foam support 05-900800	Dimensions	Height	55 mm
		Length	330 mm
		Depth	140 mm
		Hole diameter	35 mm
		Place number	21
	Material	Moltoprene, dark grey	
Foam support 05-900850	Dimensions	Height	55 mm
		Length	330 mm
		Depth	140 mm
		Hole diameter	70 mm
		Place number	5
	Material	Moltoprene, dark grey	

Date of emission: 30/08/2024

Rev. 01