

Base for Archives



In the photo, from the top:

- Cover
- Module Histoslides
- Module Histoblocks
- Module Colorteca
- Module for macrosections
- Base

CODE	DESCRIPTION	PACKAGING
03-5000-BA	Base for modular archives Histoslides, Histoblocks and Colorteca	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 22/07/2024

Rev. 04

Base for modular archives, realized in sheet metal, thickness 9/10, japanned white.
It has to be used with:

- Histoslide, code 03-5000-14,
- Histoblock, code 03-B900,
- Macrosection Istoteca, code 03-7000,
- Colorteca, series code 03-COLOR13.

Technical Specifications

Technical data	Material	Metal sheet
	Dimensions (l x p x h)	490 x 490 x 60 mm
	Electrical wires	Not necessities
Packaging	Primary container	Cardboard box
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/07/2024

Rev. 04