

Bio-Block

Archive for paraffin blocks



CODE	DESCRIPTION	PACKAGING
03-3000	Bio-Block, plastic archive for paraffin blocks	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 10/11/2022
Rev. 03

Modular archive with 8 drawers for the paraffin blocks' storage; each drawer has seven partitions.

The total Bio-Block capacity is 1.848 histology cassettes (considering 3-4 mm of paraffin per each embedding cassette) or 1.000 embedding rings.

Each module has a guide rail system all over its surface, so it is possible to arrange it in different shapes.

The Bio-Block is made in resistant plastic so it is light and handy for the shifting.

Electrical connections are not necessary.

Technical Specifications

Technical data	Material	Polystyrene
	Dimensions (l x p x h)	24 x 30 x 40 cm
	Weight	6 kg
	Drawers per module	8
	Block capacity per module	About 1.848
Packaging	Primary container	Cardboard box
	Lable	Paper adhesiv lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Plastic waste, dispose in conformity with local laws.

Date of issue: 10/11/2022

Rev. 03