



BIOPSY CASSETTES

Embedding cassettes for biopsies



CODE	COLOUR	PACKAGING	UDI-DI
07-7200	White	1.500 pcs (3 dispensers with 500 pcs each)	08033976235338
07-7210	Orange	1.500 pcs (3 dispensers with 500 pcs each)	08033976235345
07-7220	Blue	1.500 pcs (3 dispensers with 500 pcs each)	08033976235352
07-7230	Yellow	1.500 pcs (3 dispensers with 500 pcs each)	08033976235369
07-7250	Pink	1.500 pcs (3 dispensers with 500 pcs each)	08033976235376
07-7260	Green	1.500 pcs (3 dispensers with 500 pcs each)	08033976235383
07-7280	Grey	1.500 pcs (3 dispensers with 500 pcs each)	08033976235390



In vitro Diagnostic – Medical Device
IVD Class A, Reg. UE 2017/746



Manufacturer: Bio-Optica Milano S.p.A.

Basic UDI: 080339762W05030501YL



Disposable

Date of issue: 16/05/2022
Rev. 001

Cassettes for histological embedding of biopsies and small samples.

TECHNICAL FEATURES

Cassettes with attached lid and square holes with a double closure system to guarantee a perfect lock-up during the processing steps.

Made in acetalic polymer: this polymer is stable to high temperatures, microwave operations, solvents, acids and alkalis; therefore, this material is optimal for the processing and the embedding passages in the histological process.

The sample's identification is possible thanks to the wizen surface on three cassette's sides. These surfaces can be written with pencils, histology pens (11-50 or 11-400) and they're compatible with automatic cassettes printers Leica IP and Sakura Autowrite. Biopsy cassettes are manufactured with our best technical skills, therefore, due to the complexity of the details themselves and the type of material used, is to be considered a physiological quantity of defective cassettes equal to 5,000 PPM (Parts Per Million), corresponding to a percentage of 0.5%. Within this amount, Bio-Optica does not accept complaints and / or requests for warranty replacement.

Technical data

Technical specifications	Chemical composition	Acetalic polymer (polyoximethylene polyacetal).
	Internal dimensions	30 x 25 x 5 mm
	External dimensions	40 x 28 x 6 mm
	Holes dimensions	(0.6 ± 0.2) x (0.6 ± 0.2) mm
	Surface inclination	45°
	Weight	3,40 g
Packaging	Primary container	Cardboard dispensers with frontal opening to ease the cassettes extraction.
	Secondary container	Carton box.
Storage	Storing conditions	Particular conditions are not necessary for this kind of material.
	Stability	After the first opening, the product is stable.
	Validity	Not applicable to this kind of product
Warning and precautions	Instruction of use	Not requested for this product.
	Classification of the product	The product is intended for professional laboratory use for healthcare professionals. The product is not chemically dangerous. No warning measures and precautions are required.
	Disposal	Observe all state and local environmental regulations regarding waste disposal.
	Recommendations	In the event of a serious accident, we recommend that you immediately inform Bio-Optica Milano S.p.A. and the competent authorities.

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022

Date of issue: 16/05/2022

Rev. 001