

Colorteca

Archive for microscope slides and/or blocks



In the photo, from the top:

- Cover
- Module for macrosections
- Module Colorteca
- Module Istobloc
- Module Istoglass
- Base

CODE	DESCRIPTION	PACKAGING
03-COLOR13	White carrying structure	1 pcs
03-CA7100	White drawer	1 pcs
03-CA7110	Orange drawer	1 pcs
03-CA7120	Blue drawer	1 pcs
03-CA7130	Yellow drawer	1 pcs
03-CA7140	Lilac drawer	1 pcs
03-CA7150	Pink drawer	1 pcs
03-CA7160	Green drawer	1 pcs
03-CA7180	Grey drawer	1 pcs
03-CA7100S	White drawer with internal separators for 156 slides	1 pcs
03-CA7110S	Orange drawer with internal separators for 156 slides	1 pcs
03-CA7120S	Blue drawer with internal separators for 156 slides	1 pcs
03-CA7130S	Yellow drawer with internal separators for 156 slides	1 pcs
03-CA7140S	Lilac drawer with internal separators for 156 slides	1 pcs
03-CA7150S	Pink drawer with internal separators for 156 slides	1 pcs
03-CA7160S	Green drawer with internal separators for 156 slides	1 pcs
03-CA7180S	Grey drawer with internal separators for 156 slides	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 16/09/2024

Rev. 04

Colorteca is a modular archive for embedded samples and slides. Each Colorteca has 13 drawers, usable for the storing of both blocks and slides.

The plastic drawers are available in different colours to ease the blocks and slides' storage.

Colorteca has been designed for the need to have a unique archive to store and, at the same time, classify embedded blocks and slides.

A white japanned metal sheet carrying structure, thickness 10/10, with 13 coloured plastic drawers, sliding on a guide composes a Colorteca.



The drawers are available in eight different colours, with the same bio cassettes tonalities. It is possible to choose all thirteen drawers of the same colour or drawers of different colours.

Each drawer contains two separating parts that help to maintain the uniformity and the vertical slides and cassettes position.

Colorteca may be used with the base (code 03-5000-BA) and the cover (code 03-5000-CO).

Colorteca may be also overlapped to Istoglass (code 03-5000-14).

Optional accessories

CODE	DESCRIPTION	PACKAGING
03-5000-BA	Base for Istoglass, white colour	1 pcs
03-5000-CO	Cover for Istoglass, white colour	1 pcs

Technical Specifications

Technical data	Material of carrying structure	Metal sheet
	Drawers material	Plastic
	Dimensions (l x p x h)	490 x 490 x 120 mm
	Weight (empty)	6 kg
	Weight (full with slides)	32 kg
	Weight (full with cassettes)	11 kg
	Drawers per module	13
	Capacity per module (cassettes)	624 (about 48 per drawer)
	Capacity per module (rings)	312 (about 24 per drawer)
	Capacity per module (slides)	4290 (about 330 per drawer)

Date of issue: 16/09/2024

Rev. 04

Packaging	Primary container	Cardboard box
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 16/09/2024

Rev. 04