

# **HISTOCOLD**

*Storage system for a safe transport of histological specimens  
From the surgery room to the PA lab  
(In conformity with the European Regulation 605/2014)*



CODE	DESCRIPTION	PACKAGING
05-100300	Fridge for transport (65L)	1 pcs
05-100400	Cart	1 pcs
65-SL2000	UPS for fridge	1 pcs

Date of issue: 23/10/2018  
Rev. 01

### Fridge (65L)



Dimensions	l	700 mm
	d	440 mm
	h	470 mm
	Capacity	65 L
	Weight	24 kg
Characteristics	Temperature	Regolabile da +10°C a -20°C (oscillazione max +/- 1°C)
	Power supply	Direct current (battery) 12V/24V - 110V/240V / 50Hz-60Hz AC – 2 different cables according to the desired type of supply.
	Electrical power	45 W

### Cart for transport



Material	Stainless steel	
Dimensions	l	101 mm
	d	550 mm
	h	945 mm
	Weight	11 kg

Date of issue: 23/10/2018

Rev. 01