




Product Features



Precise line scanning technology ensures perfect quality at high speed



Simple operation
Automatically scans after closing the lid




Access all images via Internet on any device, phone or Ipad

Application Circumstances



Clinical Diagnosis

- Information management
- AI-assisted diagnosis



Freezing In Medical Operation

- Fast Scanning
- Timely diagnosis



Training

- Resources sharing
- Quality assurance






Remote Group Diagnosis

- Quick transmission
- Real-time image viewing

Performance Index

Acceptable Slide Size	1 inch×3 inch
Capacity	1-20 slides
Lighting	LED
Scanning Manner	line scanning
Objective Lens	OLYMPUS super apochromatic objective lens Numerical Aperture 0.8
Movement Control	magnetic linear motor control
Magnification	20 / 40
Resolution	0.25μm/pixel (20X) 0.125μm/pixel (40X)
Scanning Speed	≤25s (20x, 15mm*15mm) ≤40s (40x, 15mm*15mm)
Z-Stack Scan	yes
Focusing	automatic/manual
Precise Identification	yes
Label Identification	Bar code/QR
Image Format	KFB/TIF/DICOM
Work Station	CPU: Intel Core I7 Hard Disk: 1TB Memory: 64G Monitor: 30 inch
Dimension (L*W*H)	486*461*561mm

Product Dimension



*Unit: mm

Contact Us

Tel/WhatsApp/Wechat: +86 186-7233-4799

Email: info@kfbio.cn

Web: kfbio pathology.com

All the information in this material is for reference only, the company reserves all the rights, including updating and revision, sales representatives shall be consulted on issues of parameter updating.

