

Diamond scriber

General features

Diamond scriber for the classification of histological and cytological slides.



Code	Packaging	Length
08-DS2F	1 pcs	15 cm

Distributor: Bio-Optica Milano S.p.A.

Technical features

High quality, precision scribing pens for marking, engraving glass. To use for the classification of histological and cytological samples.

With very fine straight diamond tip (tip diameter: 1,25 mm) for the best and clear writing.

With aluminium hex handle.

Considering the product's characteristics, particular care for storing and use is not requested.

Expiry date: not applicable.