



MICROTOME BRUSHES



CODE	DESCRIPTION	PACKAGING	UDI-DI
08-0822	Small, for biopsies sections collection	4 pcs.	08033976236267
08-0823	Small, for sections collection	4 pcs.	08033976236274
08-0824	Medium, for sections collection	4 pcs.	08033976236281
08-0825	Medium, for microtome blade cleaning	2 pcs.	08033976236298
08-0826	Big, for microtome cleaning	2 pcs.	08033976236304
08-0827	Big, for cryostat cleaning	2 pcs.	08033976236311
08-0828	Set of 5 brushes for cryostat (include codes 08-0822, 08-0823, 08-0824, 08-0824 and 08-0826)	1 pc.	08033976236328
08-0829	Set of 5 brushes for microtome (include codes 08-0822, 08-0823, 08-0824, 08-0824 and 08-0825)	1 pc.	08033976236335



In vitro Diagnostic – Medical Device
IVD Class A, Reg. UE 2017/746



Manufacturer: Bio-Optica Milano S.p.A.

Basic UDI: 080339762W05030580ZA

Date of issue: 16/05/2022
Rev. 001

Brushes to use in microtomy and cryo-microtomy.

Technical data

	Code	Tip	Handle
Material	08-0822, 08-0823, 08-0824	Cow fur	Wood
	08-0825, 08-0826	Bristle (30% animal 70% synthetic)	Polypropylene
	08-0827	Bristle (30% animal 70% synthetic)	Wood
Packaging	Primary container	Plastic bag	
	Label	Paper label	
	Secondary container	Not required	
Storage	Storing condition	None	
	Storage temperature	There are no particular limitations.	
	Cleaning method	Clean with cloth or paper. Don't use solvents.	
	Stability	After the first opening, the product is stable.	
	Validity	Not applicable.	
Warning and precautions	Instruction for use	Not provided for this product.	
	Classification of the product	The product is intended for professional laboratory use for healthcare professionals. The product is not chemically dangerous. No warning measures and precautions are required.	
	Disposal	Observe all state and local environmental regulations regarding waste disposal.	
	Recommendations	In the event of a serious accident, we recommend that you immediately inform Bio-Optica Milano S.p.A. and the competent authorities.	

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022

Date of issue: 16/05/2022
Rev. 001