



## TRANSPORT BAGS



CODE	DESCRIPTION	PACKAGING	UDI-DI
44-9590	Transport Bags	500 pcs	08033976239534



In vitro Diagnostic – Medical Device  
IVD Class A, Reg. UE 2017/746



Manufacturer: Bio-Optica Milano S.p.A.

Basic UDI: 080339762W05020101XM



Disposable

Date of issue: 16/05/2022  
Rev. 001

Bags for transport of biological samples (anatomical or cytological liquids).

### GENERAL FEATURES

The bags for the transport of samples potentially infected have a closure with pressure zip and an external pocket for documents that allows protecting patient's documentation in case of container's spills.  
The bags have on external side Biohazard symbol and the use's instructions.

### Technical data

Technical data	Material	LDPE (Low Density Polyethylene)
	Total bag dimensions	160 x 270 mm
	Useful bag dimensions (until zip)	160 x 240 mm
	Pocket dimensions	160 x 180 mm
	Tolerate temperatures	-10 / +60°C
Packaging	Primary container	Carton box
	Label	Paper label
	Secondary container	Not requested
Storage	Storing condition	None
	Storage temperature	There are no particular limitations.
	Stability	After the first opening, the product is stable.
	Validity	Not applicable.
Warning and precaution	Classification of the product	The product is intended for professional laboratory use for healthcare professionals. The product is not chemically dangerous. No warning measures and precautions are required.
	Disposal	Observe all state and local environmental regulations regarding waste disposal.
	Recommendations	In the event of a serious accident, we recommend that you immediately inform Bio-Optica Milano S.p.A. and the competent authorities.

REVISION N°	REASON	REVISION DATE
001	Regulation adjustment UE 2017/746 - IVDR	16/05/2022

Date of issue: 16/05/2022  
Rev. 001