



BioCassettes





Cassettes with attached lid and rectangular holes, with a double closure system to guarantee a perfect lock-up during the processing phase.

The sample's identification is possible thanks to the wizen surface on three cassette's sides. These surfaces can be written with pencils, histology pens and they're compatible with automatic cassettes printers Leica IP and Sakura Autowrite.

Bio Cassettes

Acetal resin cassettes for embedding standard samples with lid or with separate lid. Available in 8 different color:

- White 07-7100 or 07-8100 (separate lid)
- Orange 07-7110 or 07-8110 (separate lid)
- Blue 07-7120 or 07-8120 (separate lid)
- Yellow 07-7130 or 07-8130 (separate lid)
- Lilac 07-7140 or 07-8140 (separate lid)
- Pink 07-7150 or 07-8150 (separate lid)
- Green 07-71060 or 07-8160 (separate lid)
- Grey 07-7180



Biopsy Cassettes

Acetal resin cassettes for embedding biopses samples with lid or with separate lid. Available in 8 different color:

- White 07-7200 or 07-8200 (separate lid)
- Orange 07-7210 or 07-8210 (separate lid)
- Blue 07-7220 or 07-8220 (separate lid)
- Yellow 07-7230 or 07-8230 (separate lid)
- Lilac 07-7240 or 07-8240 (separate lid)
- Pink 07-7250 or 07-8250 (separate lid)
- Green 07-7260 or 07-8260 (separate lid)
- Grey 07-7280 or 07-8280 (separate lid)



Mega Cassettes

Double height cassettes with lid. Available in white, code 07-7300.



Super Mega Cassettes

Cassettes for large samples with separate lid, designed with larger mesh to increase the contact surface of the paraffin and prevent the sample from detaching during microtome cutting.

Available in white, code 07-7000.

Stacked Bio Cassettes with lid with tape for printers

Standard Embedding cassettes with double closure system to guarantee a perfect resistance during tissue processing.

Suitable surface for automatic cassette printers. The cassettes are stacked with tape, with attached lid.

Available in 7 different color:

- White 07-9700
- Orange 07-9710
- Blue 07-9720
- Yellow 07-9730
- Lilac 07-9740
- Pink 07-9750
- Green 07-9760

Compatible with these writing system:

- Primera
- Leica IPC
- Sakura Smart Write
- Sakura Auto Write