


## Disposable Microtome Blades

- › S35
- › R35
- › N35
- › N35HR
- › A35
- › A22
- › S22
- › C35
- › S35L
- › S35LL
- › High Profile
- › SH35W



pfm medical ag offers a broad range of FEATHER® disposable microtome blades for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

[www.pfmmedical.com](http://www.pfmmedical.com)



pfm medical ag offers a broad range of **FEATHER® disposable microtome blades** for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

**S35**



**Use**

The S35 disposable microtome blade is a standard blade for routine practice and biopsies.

**Details**

- ▶ With excellent sectioning quality from the beginning
- ▶ Very long service life
- ▶ Short trimming time

**Technical Data**


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type S35**

**REF 207500000**

**R35**



**Use**

The R35 disposable microtome blade is extra suitable for routine ribbon sectioning and for hard tissues.

**Details**

- ▶ Specially hardened cutting edge
- ▶ Short trimming time

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type R35**

**REF 207500005**

## N35



### Use

The N35 disposable microtome blade is particularly suitable for hard tissues.

### Details

- ▶ Longer service life
- ▶ Suitable for all routine sectioning, also for soft tissue and biopsies
- ▶ Also suitable for ribbon sectioning after trimming
- ▶ Low compression

### Technical Data


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

#### Type N35

REF 207500006

## N35HR



### Use

The N35HR disposable microtome blade is particularly suitable for sectioning hard and difficult tissues.

### Details

- ▶ Suitable for routine and ribbon sectioning
- ▶ Low compression
- ▶ Longer service life
- ▶ Short trimming time
- ▶ New coating
- ▶ Specially finished cutting edge

### Technical Data


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

#### Type N35HR

REF 207500014

**A35**



**Use**

The A35 disposable microtome blade allows the production of very thin sections and is suitable for hard and soft tissues.

**Details**

- ▶ Very long service life (also for hard tissue)
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type A35**

**REF 207500011**

**A22**



**Use**

The A22 disposable microtome blade allows the production of very thin sections in routine use and research and is suitable for fibrous tissues (biopsies/fabrics).

**Details**

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp
- ▶ Also suitable for ribbon sectioning after trimming

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type A22**

**REF 207500002**

**S22**



**Use**

The S22 disposable microtome blade allows the production of extra-thin sections and is suitable for soft tissues/fabrics.

**Details**

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type S22**

**REF 207500001**

**C35**



**Use**

The C35 disposable microtome blade is specifically designed for cryostat sectioning.

**Details**

- ▶ Carbon steel blade
- ▶ Short trimming time

**Technical Data**

- ▶ **Material:** carbon steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

**Ordering Information**

**Type C35**

**REF 207500003**

## S35L



### Use

The S35L disposable microtome blade with a length of 120 mm is particularly suitable for large tissue blocks and routine sectioning.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 120 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

### Ordering Information

**Type S35L**

**REF** 207500004

## S35LL



### Use

The extra-long S35LL disposable microtome blade with a length of 180 mm is suitable for large tissue blocks and routine sectioning.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 180 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

### Ordering Information

**Type S35LL**

**REF** 207500009

## High Profile



### Use

The High Profile disposable microtome blade is suitable for all paraffin embedded tissue blocks in routine applications.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 14 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

**Type High Profile**

**REF 207500010**

## SH35W



### Use

The SH35W disposable microtome blade is particularly suited for sectioning resin embedded tissue and hard materials. Suitable for use with the pfm medical disposable blade holder HP TC (REF 030065).

### Details

- ▶ High Profile blade with Tungsten Carbide edge

### Technical Data

- ▶ **Material:** stainless steel with Tungsten Carbide edge
- ▶ **Length:** 80 mm
- ▶ **Height:** 15 mm
- ▶ **PU:** 3 pcs.

### Ordering Information

**Type SH35W**

**REF 207500015**

## Blade Depot



### Use

The blade depot provides suitable storage for FEATHER® disposable microtome blades.

### Details

- ▶ Storage capacity: approx. 200 blades
- ▶ Accessorie: protective cover blade deposit

### Manufacturer

pfm medical ag  
Wankelstraße 60  
50996 Köln, Germany

### Ordering Information

#### Blade Depot

REF 030000

#### Protective cover blade deposit

REF 030100

### Ordering Information

#### Disposable Microtome Blades

| Type         | Colour coding | REF       | Material                                   | Length | Height | Blade angle | PU      |
|--------------|---------------|-----------|--|--------|--------|-------------|---------|
| S35          | dark blue     | 207500000 | stainless steel                            | 80 mm  | 8 mm   | 35°         | 50 pcs. |
| R35          | black         | 207500005 | stainless steel                            | 80 mm  | 8 mm   | 35°         | 50 pcs. |
| N35          | light green   | 207500006 | stainless steel                            | 80 mm  | 8 mm   | 35°         | 50 pcs. |
| N35HR        | dark green    | 207500014 | stainless steel                            | 80 mm  | 8 mm   | 35°         | 50 pcs. |
| A35          | light blue    | 207500011 | stainless steel                            | 80 mm  | 8 mm   | 35°         | 50 pcs. |
| A22          | pink          | 207500002 | stainless steel                            | 80 mm  | 8 mm   | 22°         | 50 pcs. |
| S22          | grey          | 207500001 | stainless steel                            | 80 mm  | 8 mm   | 22°         | 50 pcs. |
| C35          | white         | 207500003 | carbon steel                               | 80 mm  | 8 mm   | 35°         | 20 pcs. |
| S35L         | blue          | 207500004 | stainless steel                            | 120 mm | 8 mm   | 35°         | 20 pcs. |
| S35LL        | green         | 207500009 | stainless steel                            | 180 mm | 8 mm   | 35°         | 20 pcs. |
| High Profile | brown         | 207500010 | stainless steel                            | 80 mm  | 14 mm  | 35°         | 50 pcs. |
| SH35W        | blue          | 207500015 | stainless steel with Tungsten Carbide edge | 80 mm  | 15 mm  | 35°         | 3 pcs.  |

### Contact

Should you have any questions our Customer Solutions Team is happy to advise you.

✉ service@pfmmedical.com

☎ +49 (0)2236 9641-220

☎ +49 (0)2236 9641-51

pfm medical ag  
Wankelstraße 60  
50996 Köln, Germany

Exclusive Distributor  
in Europe

Certified according to  
DIN EN ISO 13485

#### Manufacturer

FEATHER® Safety Razor Co. Ltd. Registered Trademark in Japan, Products are subject to alterations,

