




Biocare Medical LLC

EC Declaration of Conformity

Date of Issue:	February 3, 2021
Certificate Reference:	All Biocare IVD products are self-declared for compliance to Annex III of the IVDD. None contain any analytes listed in Annex II List A or B
Directive:	98/79/EC IVD Directive of 27 October 1998
Conforming Product:	See attached Product List
Manufacturer:	Biocare Medical LLC 60 Berry Drive, Pacheco, CA 94553, U.S.A
Authorized Representative:	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands +31.70.345.8570 - phone +31.70.346.7299 - fax europe@emergogroup.com

Biocare Medical declares that the above-mentioned products meet the provision of the Council Directive 98/79/EC for In Vitro Diagnostic Medical Devices and Directive 98/79/EC as transposed in the national laws of the Member States.

Signed: 
Johanna Saito, Regulatory Affairs Manager




Biocare Medical LLC

Dichiarazione di conformità CE

Traduzione in italiano

Data di rilascio:	3 febbraio 2021
Riferimento del certificato:	Tutti i prodotti Biocare IVD sono auto-dichiarati per la conformità all'allegato III dell'IVDD. Nessuno contiene analetti elencati nell'allegato II, elenco A o B
Direttiva:	Direttiva IVD 98/79/CE del 27 ottobre 1998
Prodotto conforme:	Vedi l'elenco dei prodotti allegato
Fabbricante:	Biocare Medical LLC 60 Berry Drive, Pacheco, CA 94553, USA
Rappresentante autorizzato:	Emergo Europe Prinsessegracht 20 2514 AP L'Aia Paesi Bassi +31.70.345.8570 - telefono +31.70.346.7299 - fax europe@emergogroup.com

Biocare Medical afferma che i suddetti prodotti soddisfano le disposizioni della direttiva 98/79/CE del Consiglio per i dispositivi medico-diagnostici in vitro e della direttiva 98/79/CE recepita nelle legislazioni nazionali degli Stati membri.

Firmato: 
Johanna Saito, Responsabile Affari Regolatori



Biocare Medical LLC

Part Number	Description
TWN20	20% Tween 20
GM601	4plus Biotinylated Goat Anti-Mouse IgG
GR602	4plus Biotinylated Goat Anti-Rabbit IgG
GR608	4plus Biotinylated Goat Anti-Rabbit IgG (Mouse and rat adsorbed)
GU600	4plus Biotinylated Universal Goat Link
HP504UM	4plus HRP Universal Detection
HP504US	4plus HRP Universal Detection
AP605	4plus Streptavidin AP Label
HP604	4plus Streptavidin HRP Label
AP506US	4plus Universal AP Detection
VLTR3138	Adipophilin
ACI/API3138	Adipophilin
VLTM3041	ALK [5A4]
ACI/API3041	ALK [5A4]
CP/PP028	Alpha-1-Fetoprotein (AFP)
CM/PM125	Amyloid A
PP132	Amyloid P
ACI/API109	Androgen Receptor
CB910	Antigen Decloaker, 10X
OPRI6006	Antigen Retrieval 1 (AR1)
ORI6006	Antigen Retrieval 1 (AR1), high pH
ORI6005	Antigen Retrieval 2 (AR2), low pH
ADP1002	Aqua DePar, 10X
VLTR3058	Arginase-1
ACI/API/ OAI / AVI3058	Arginase-1
TWA20	Automation Tween 20, 20X
AB972	Avidin-Biotin Kit
BE965	Background Eraser
BP974	Background Punisher
BS966	Background Sniper
BT967	Background Terminator
BJP811	Bajoran Purple Chromogen Kit
ACI/API/VLTR3247	BAP1
ACI/API/VLTM3241	BCA-225 [Cu-18]
OPAI003	Bcl-2
AVI003	Bcl-2
VLTM003	Bcl-2



Biocare Medical LLC

CM/PM/IP/OAI003	Bcl-2
AVI410	Bcl-6 [LN22]
VLTM410	Bcl-6 [LN22]
CM/PM/OAI410	Bcl-6 [LN22]
OPAI410	Bcl-6 [LN22]
AVI107	Ber-EP4
PM/OAI107	Ber-EP4
API3112	Ber-EP4 + BG8
DS900	Betazoid DAB Buffer
BDB2004	Betazoid DAB Chromogen Kit
CM/PM418	BOB-1
BD1000	Borg Decloaker, RTU
CM345	BRCA-1
AVI3203DSK	Breast Cocktail (CK HMW/p63 + CK7/8/18)
ACI/API3134	C4d (RM)
CM/PM101	CA 125
AVI123	CA 19-9
CM/PM123	CA 19-9
VLTM101	CA125
CP/PP072	Calcitonin
VLTR072	Calcitonin
VLTM172	Calponin
CM/PM172	Calponin
OPAI092	Calretinin
CP/PP/IP/OAI092	Calretinin
VLTR092	Calretinin (Rb)
CM/PM058	Carcinoembryonic Antigen (CEA {M})
CP/PP/IP009	Carcinoembryonic Antigen (CEA {P})
DBC859	Cardassian DAB Chromogen Kit
TRP955	Carezyme I: Trypsin Kit
PEP956	Carezyme II: Pepsin
PRT957	Carezyme III: Pronase Kit
CP/PP229	Caspase-3 (Cleaved)
CATHE	CAT Hematoxylin
ACI/API144	Cat Scratch (Bartonella henselae)
VLTM129	CD10
OPAI129	CD10
CM/PM/IP/OAI129	CD10
API3212	CD10 [UMAB236]



Biocare Medical LLC

ACI/API3117	CD103 (RM)
CME/PME/IP/OAI296	CD117/c-kit
OPAI296	CD117/c-kit
VLTR296	CD117/c-kit (Rb)
ACI/API3122	CD11c (Leu-M5)
ACI/API3201	CD137
VLTM167	CD138
CM/PM/IP167	CD138
API3193DS	CD138 + Cyclin D1
API3169DS	CD138 + Ki-67
CM/PM029	CD15 [MMA]
CM/PM/IP/OAI073	CD15 Cocktail
CM/PM/OAI353	CD163
CM/PM310	CD19
VLTM3158	CD1a
ACI/API3158	CD1a [O10]
OPAI004	CD20 [L26]
VLTM004	CD20 [L26]
CM/PM/IP/OAI004	CD20 [L26]
CM/PM/OAI142	CD21
CM/PM169	CD22
ACI/API3229	CD22
VLTM100	CD23
CM/PM/OAI100	CD23
CME/PME324	CD3
OPAI324	CD3
VLTM110	CD3 (M)
VLTR215	CD3 (P) (RP)
ACI/API3170	CD3 [BC33]
OPAI3170	CD3 [BC33]
VLTM3152	CD3 [LN10]
ACI/API3152	CD3 [LN10]
CM/PM 110	CD3 T-Cell (M)
OPAI110	CD3 T-cell (M)
CP/PP215	CD3(P)
PM/IP031	CD30 (Ki-1)
CM/PM346	CD30 (Ki-1)
VLTM074	CD30 Cocktail
PM074	CD30 Cocktail



Biocare Medical LLC

CM/PM/OAI131	CD31
CM/PM347	CD31 (PECAM-1)
VLTM3116	CD33
ACI/API3116	CD33
OPAI084	CD34
VLTM084	CD34
CM/PM/IP/OAI084	CD34
AVI167	CD38
CM/PM153	CD4 (Helper/Inducer)
VLTR3209	CD4 (RM)
API3209	CD4 (RM)
VLTM3148	CD4 [4B12]
ACI/API3148	CD4 [4B12]
API3157DSAA	CD4 + CD8
CM/PM/IP005	CD43
CM/PM318	CD44
PM380	CD44
CM/PM006	CD45RO [UCHL-1]
CM/PM/OAI099	CD5 (M)
CM/PM/OAI164	CD56
CM/PM/OAI007	CD57 (Natural Killer Cell)
AVI033	CD68 [KP1]
OPAI033	CD68 [KP1]
CM/PM/IP/OAI033	CD68 [KP1]
VLTM158	CD7
CM/PM/OAI158	CD7
ACI/API3110	CD71
VLTM067	CD79a
CM/PM067	CD79a
VLTR311	CD8
CRM/PRM311	CD8
ACI/API3160	CD8 [C8/144B]
OPAI3160	CD8 [C8/144B]
ACI/API3219	CD8 [C8/1779R]
VLTM3160	CD8 [CD8/144B]
VLTR392	CD99
OPAI392	CD99
CME/PME/ OAI392	CD99
VLTM3111	CDH17



Biocare Medical LLC

ACI/API3111	CDH17 (M)
OPAI226	CDX2
CM/PM/IP/OAI/ VP226	CDX2
API3135DS	CDX2 (M) + CDH17 (RM)
ACI/API3144	CDX2 (RM)
VLTR3144	CDX2 (RM)
ACI/API AVI3184	CDX2 [BC39]
PM367DS	CDX2 + CK7
VLTM058	CEA (M)
ORI6031K	Chamber Cleaning Kit
VLTM010	Chromogranin A
AVI010	Chromogranin A
OPAI010	Chromogranin A
CM/PM/IP OAI010	Chromogranin A
VLTR339	CK 7
API/IPI/ OAI3154DS	CK HMW + p63 + AMACR (RM)
OPAI3154	CK HMW + p63 + AMACR (RM)
VLTMR3154	CK HMW + p63 + AMACR (RM)
OAI3124K	CK HMW + p63, 2X
VLTM056	CK LMW 8/18
VLTM3061	CK18
VLTM062	CK20
OPAI235	CK5 + p63
PM235	CK5 + p63
PM/OAI360DS	CK5/14 + p63 + CK7/18
VLTMR360	CK5/14 + p63 + CK7/18
OPAI360	CK5/14 + p63 + CK7/18
VLTM061	CK7
CRM/PRM/IP339	CK7
API3161	CK8/18 (RM)
ACI/API3121	Claudin-4
ACI218	Clusterin
AVI415	c-Myc
CME/PME415	c-Myc
CM/PM112	Collagen IV
CRM/PRM306	COX-2



Biocare Medical LLC

ACI/API3211	CTLA-4
ACI/API3220	CTLA-4 [CAL49]
CRM/PRM/IPI307	Cyclin D1
CME/PME/ OAI432	Cyclin D1
OPAI432	Cyclin D1
OPAI307	Cyclin D1
VLTR432	Cyclin D1 (Rb)
PM081	Cytokeratin [AE1] LMW
CM185	Cytokeratin 14 (CK14)
PM176	Cytokeratin 17 (CK17)
ACI/API3061	Cytokeratin 18 (CK18)
CM/PM/OAI242	Cytokeratin 19 (CK19)
OPAI062	Cytokeratin 20 (CK20)
CM/PM/IP/OAI062	Cytokeratin 20 (CK20)
CM/PM/OAI234	Cytokeratin 5 (CK5)
CME/PME430	Cytokeratin 5 (CK5)
ACI/API/OAI3025	Cytokeratin 5/14 Cocktail
OPAI105	Cytokeratin 5/6 (CK 5/6)
AVI105	Cytokeratin 5/6 (CK 5/6)
CM/PM/IPI105	Cytokeratin 5/6 (CK 5/6)
OPAI061	Cytokeratin 7 (CK7)
CM/PM/IP/OAI061	Cytokeratin 7 (CK7)
VLTM127	Cytokeratin HMW
CM/PM/IP/OAI127	Cytokeratin HMW (34βE12)
AVI127	Cytokeratin HMW [34βE12]
OPAI127	Cytokeratin HMW [34βE12]
OPAI056	Cytokeratin LMW (8/18)
CM/PM/IP/OAI056	Cytokeratin LMW 8/18 (LMW CK 8/18)
OPAI266	D2-40 (Lymphatic Marker)
CM/PM/IP/OAI266	D2-40 (Lymphatic Marker)
PD900	Da Vinci Green Diluent
DA000-250-Kit	DAB Away
OPRI6055	DAB Buffer
OPRI6054	DAB Chromogen
DB801	DAB Chromogen Kit
DS830	DAB Sparkle
DS854	DAB Substrate Buffer
DC2012-220V	DC2012 Decloaking Chamber NxGen (220v model)
DC2008INTL	Decloaking Chamber™ Plus



Biocare Medical LLC

DC2002INTL	Decloaking Chamber™ Pro
BRI4015	Deep Space Black™ Chromogen Kit
DNS001	Denaturing Solution (Elution Step)
BRI4044	DepART Solution
DRY2008INTL	Desert Chamber Pro
CM/PM/IP/OAI036	Desmin
CM/PM419	Desmoglein 3
ACI/API3018	Desmoglein 3 + CK5
PPM428DS	Desmoglein 3 + Napsin A
API3067	Desmoglein 3 + p40 (M)
API3132DS	Desmoglein 3 + p40 (M) + Napsin A (RM)
OPRI6001	Dewax Solution 1 (DS1)
OPRI6003	Dewax Solution 2 (DS2)
OPRI6004K	Dewax Solution Kit
ORI6004K	Dewax Solution Kit
BRI2006	Diva Decloaker, 100X
DV2004X	Diva Decloaker, 10X
DV2005	Diva Decloaker, 20X
DV2004	Diva Decloaker, RTU
BRI4036	DNA Hybridization Buffer
CM/PM/OAI385	DOG1
ORI6049	DS Enzyme
OPRI6049	DS Enzyme
OPRI6060K	DS Enzyme
ACI/API3230	E-Cadherin
CM/PM/IP170	E-Cadherin
ACI/API3012	E-Cadherin (RM)
EM897	EcoMount
HTE	Edgar Degas Eosin
CB917	EDTA Decloaker, 5X
ACI063	Epidermal Growth Factor Receptor (EGFR)
CM/PM143	Epithelial Membrane Antigen (EMA [Mc-5])
ACI/API3038	Epithelial Membrane Antigen (EMA) [E29]
API111	Epstein-Barr Virus (EBV)
ACI3147	ERCC1
CM/PM/OAI421	ERG
OPAI421	ERG



Biocare Medical LLC

API3017	ERG (M), 2X
API3176DS	ERG + C4d
APR3176DS	ERG + C4d
API437DS	ERG-2 (ERG + CK5)
VLTMR437	ERG-2 (ERG + CK5)
ACI/API054	Estrogen Receptor (ER) [1D5]
API308	Estrogen Receptor (ER) [6F11 + SP1]
ACI/API093	Estrogen Receptor (ER) [6F11]
ACI/API/OAI301	Estrogen Receptor (ER) [SP1]
CP/PP039	Factor VIII
CM/PM/IP/VP357	Factor XIIIa
OPAI357	Factor XIIIa
OPRI6041	Fast Red Buffer
OPRI6040	Fast Red Chromogen
ORI6042K	Fast Red Chromogen Kit
FB813	Ferangi Blue™ Chromogen Kit 2
FAD901	Fluorescence Antibody Diluent
FP001	Fluoro Care Anti-Fade Mountant
BRI4006K	Folate Receptor alpha IHC Assay Kit
CM/PM411	Follicle Stimulating Hormone (FSH)
API3164	FOXP3
ACI/API3197	FOXP3 [86D]
ACI/API3174	Galectin-3
CM/PM/OAI405	GATA-3
VLTM405	GATA-3
OPAI405	GATA-3
API/AVI3173	GATA-3 + Uroplakin II
PM317DS	GCDFP-15 + Mammaglobin
VLTM065	GFAP [GA-5] (M)
ACI/API3221	GITR [CAL8]
CP/PP040	Glial Fibrillary Acidic Protein (GFAP {P})
CM/PM065	Glial Fibrillary Acidic Protein (GFAP{M})
CM/PM408	GLUT-1
ACI/API3009	Glutamine Synthetase
VLTM396	Glypican-3
CM/PM396	Glypican-3
ACI/API3202	Granzyme B
CM/PM113	Gross Cystic Disease Fluid Protein-15 (GCDFP-15) (Brst-2)
VLTM383	H. pylori



Biocare Medical LLC

API/OAI3191	H. pylori (RM)
API3177	HBME-1
CM407	Heat Shock Protein 70
OPAI383	Helicobacter pylori
CM/PM/IP/ OAI 383	Helicobacter pylori
CM/PM/OAI166	Hepatocyte Specific Antigen (HSA)
API108	Herpes Simplex Virus 1&2 (HSV 1&2)
CME349	HIF-1 alpha
AVI057	HMB45
CM/PM/IP/OAI057	HMB45
OPAI057	HMB45
VLTM057	HMB45
CM/PM/IPI/VP165	HMB45 + MART-1 + Tyrosinase
VLTM165	HMB45 + MART-1 + Tyrosinase
OPAI165	HMB45 + MART-1 + Tyrosinase
OPAI210	HMW CK Cocktail + p63 (Basal Cell Cocktail)
CM/PM210	HMW CK Cocktail + p63 (Basal Cell Cocktail)
HTR1001	Hot Rinse, 25X
ACI/API177	HPV Cocktail Broad Spectrum (HPV-1,6,11,16,18 and 31)
ACI186	HPV-16 (CAMVIR-1)
CP/PP124	Human Chorionic Gonadotropin (hCG) beta
ACI/API3185	Human IgG
IED1203,1204	I.E.D Unit
ACI/API3210	IDO1
CM414	IGF-1R
ACI/API3115	IgG4 (M)
PWB941	Immunocare PBS, 10X
TWB943	Immunocare TBS Wash Buffer, 10X
API3180	IMP3 (RM)
CM/PM171	Inhibin, Alpha
IPS0001INTL	intelliPATH FLX Automated Slide Stainer
IPK5010	intelliPATH FLX™ DAB Chromogen Kit
IPK5017	intelliPATH FLX™ Fast Red Chromogen Kit
IPE5007	intelliPATH FLX™ Pepsin



Biocare Medical LLC

IPB5000	intelliPATH FLX™ Peroxidase Blocking Reagent
IPK5014	intelliPATH FLX™ Pronase Kit
IPK5003	intelliPATH FLX™ Universal HRP Detection Kit, 80-100 tests
IP498	intelliPATH FLX™ Universal Negative Control
IPI4006K	intelliPATH™ Folate Receptor alpha IHC
IP974	intelliPATH™ Background Punisher
IPK5027	intelliPATH™ Ferangi Blue™ Chromogen Kit
IPCS5006	intelliPATH™ Hematoxylin
IPSC5004	intelliPATH™ Multiplex Secondary Reagent 2
IPI5034	intelliPATH™ Negative Control (ER/PR)
IPK5011	intelliPATH™ Universal HRP Detection Kit
IPK5024	intelliPATH™ Warp Red Chromogen Kit
IPA5018	intelliPrep Solution
IQ1000INTL	IQ Kinetic Slide Stainer
IQ2000INTL	IQ Kinetic Slide Stainer
IQ3000INTL	IQ Kinetic Slide Stainer
IQ1000INTL-NS	IQ Kinetic Slide Stainer Without Orbital Shaker
IQ2000INTL-NS	IQ Kinetic Slide Stainer Without Orbital Shaker
IQ3000INTL-NS	IQ Kinetic Slide Stainer Without Orbital Shaker
ORI6020K	ISH Dewax Solution Kit
ORI6047K	ISH HRP Detection for ONCORE
BRI4050K	ISH HRP Detection Kit for Digoxigenin-Labeled DNA Probes
ORI6021	ISH Retrieval (AR3)
ORI6039K	ISHzyme ONCORE
VLTM3149	Kappa
VLTM3159	Kappa (M) + Lambda (P)
API3159DS	Kappa (M) + Lambda (P)
ACI/API3149	Kappa Light Chain [L1C1]
AVI325	Ki-67



Biocare Medical LLC

CRM/PRM/OAI325	Ki-67
VLTR325	Ki-67
OPAI325	Ki-67
AVI3156	Ki-67 [MIB-1]
API3156	Ki-67 [MIB-1]
VLTM3156	Ki-67 [MIB-1]
OPAI3156	Ki-67 [MIB-1]
PPM240DS	Ki-67 + Caspase-3
API3192DS	Ki-67 + CK7/8/18
API3198DS	Ki-67 + pHH3
PM237DS	L26 + CD3
ACI/API3213	LAG3 [CAL26]
VLTM3063	Lambda
ACI/API3063	Lambda Light Chain [N10/2]
VLTM016	LCA
CM/PM/IP/OAI016	Leukocyte Common Antigen (LCA) Cocktail
M1U539	MACH 1 Universal HRP-Polymer Detection Kit
MRCT523	MACH 2 Double Stain 1
MRCT525	MACH 2 Double Stain 2
MALP521	MACH 2 Mouse AP-Polymer
MHRP520	MACH 2 Mouse HRP-Polymer
RALP525	MACH 2 Rabbit AP-Polymer
RHRP520	MACH 2 Rabbit HRP-Polymer
M2U522	MACH 2 Universal HRP-Polymer Detection
M3M532	MACH 3 Mouse AP-Polymer Detection
MH540	MACH 3 Mouse HRP-Polymer
M3M530	MACH 3 Mouse HRP-Polymer Detection
M3R533	MACH 3 Rabbit AP-Polymer Detection
M3R531	MACH 3 Rabbit HRP Polymer Detection
BRI4005	MACH 4 Universal AP Detection System
M4U536	MACH 4 Universal AP Polymer Kit
M4U534	MACH 4 Universal HRP-Polymer Detection
BRI4012	MACH 4 Universal HRP-Polymer Detection System
VLTM269	Mammaglobin



Biocare Medical LLC

PM/OAI269	Mammaglobin (M)
API3186DS	MART-1 + Tyrosinase + pHH3
API3165	MART-1 + Tyrosinase + SOX10
AVI3216	MART-1 + Tyrosinase + SOX10
OPAI077	MART-1 Cocktail
CM/PM/IP/OAI077	MART-1 Cocktail
VLTM077	MART-1 Cocktail
ACI/API3131	MASH1
API3181	MCM2 + TOP2A
OPAI3181	MCM2 + TOP2A
AVI3114	Melan A (M)
ACI/API3114	Melan A (M)
VLTM3114	Melan A (M)
VLTM078	Melanoma Cocktail
CM/PM078	Melanoma Cocktail
ACI/API3175	Mesothelin
CP290	Microglia (Iba1)
VLTM423	Microphthalmia Transcription Factor (MiTF)
AVI423	Microphthalmia Transcription Factor (MiTF)
CM/PM423	Microphthalmia Transcription Factor (MiTF)
PM423H	Microphthalmia Transcription Factor (MiTF)
OPAI220	MLH-1
CM/PM/IP/OAI/AVI220	MLH-1
VLTM3214	MLH-1 [BC23]
AVI3214	MLH-1 [BC23]
ACI/API3214	MLH-1 [BC23]
AVI403	MOC-31
CM/PM403	MOC-31
PD901	Monet Blue Diluent
ORI6050	Mouse Amp HRP Detection
OPRI6050	Mouse Amp HRP Detection
BRI4049	Mouse anti-Digoxigenin
ORI6044	Mouse AP Detection
OPRI6044	Mouse AP Detection
ORI6007	Mouse HRP Detection



Biocare Medical LLC

OPRI6007	Mouse HRP Detection
VLTM219	MSH2
OPAI219	MSH2
CM/PM/OAI/AVI219	MSH2
OPAI265	MSH6
CM/PM/IP/OAI/AVI265	MSH6
VLTM3215	MSH6 [BC19]
AVI3215	MSH6 [BC19]
ACI/API3215	MSH6 [BC19]
CM/PM319	MUC-1
CM326	MUC-4
CM231	Mucin 5AC (Gastric Mucin)
OPRI6061	Multiplex Decton 1
OPRI6045	Multiplex Detection 2
ORI6045	Multiplex Detection 2
CRM/PRM/ OAI352	MUM-1
OPAI352	MUM-1
VLTR352	MUM-1
API3194DS	MUM-1 + CD56
VLTM079	Muscle Specific Actin
CM/PM/IP/OAI079	Muscle Specific Actin (MSA)
ACI/API140	Mycobacterium tuberculosis (TB)
VLTR023	Myeloperoxidase
PP023	Myeloperoxidase
VLTM115	Myogenin
CM/PM115	Myogenin
VLTM388	Napsin A
AVI434	Napsin A
CM/PM/ IPI/ OAI388	Napsin A
PP434	Napsin A
OPAI388	Napsin A
AVI388	Napsin A
ACI/API3043	Napsin A (RM)
VLTR3043	Napsin A (RM)
VLTR434	Napsin-A
ACI369	Nerve Growth Factor Receptor
CM/PM066	Neurofilament
CP/PP422	NKX3.1
OPAI422	NKX3.1



Biocare Medical LLC

ACI/API3189	NKX3.1 (RM)
CB911	Nuclear Decloaker, 10X
VLTM313	Oct 3/4
PM313	OCT-¾
ONC0001-220V	ONCORE Automated Slide Staining System (220V Model)
ORI6056K	ONCORE DAB Chromogen
OPRI6005	ONCORE PRO Antigen Retrieval 2 (AR2)
ONCPRO0001	ONCORE PRO Automated Slide Staining System
OPRI6056K	ONCORE Pro DAB Chromogen
OPRI6042K	ONCORE Pro Fast Red Chromogen
API3011DS	p120 + E-cadherin
VLTM3011	p120 + E-cadherin
ACI/API3008	p120 Catenin
VLTM3008	p120 Catenin
API3246DS	p16 + Ki-67
VLTM3246	p16 + Ki-67
API3256DS	p16 + PRAME
ACI/API/AVI3231	p16 INK4a
VLTM3231	p16 INK4a [BC42]
OPAI3231	p16 INK4a [BC42]
CM354	p21
ACI/API/AVI/IPI/ OAI3066	p40 (M)
VLTM3066	p40 (M)
API3079	p40 (M), 3X (Prostate)
ACI/API3030	p40 (P)
OPAI3030	p40 (P)
ACC3066	p40 [BC28] BSA Free
API3172	p40 + Uroplakin II
AVI298	p53
VLTR298	p53
CME/PME298	p53
OPAI298	p53
VLTM042	p53 Tumor Suppressor Protein (M)
CM/PM042	p53 Tumor Suppressor Protein (M)
VLTM163	p63



Biocare Medical LLC

CM/PM/OAI/IP/VP163	p63
VLTM391	p63 + CK5
PM391DS	p63 + CK5
VLTM427	p63 + TRIM29
PPM427DS	p63 + TRIM29
API3070	p63, 2X (Lung)
PM366	p63-2X
VLTM011	Pan CK [AE1/AE3]
OPAI011	Pan Cytokeratin [AE1/AE3]
CM/PM/IPI/VP011	Pan Cytokeratin [AE1/AE3]
CM/PM/IP/VP043	Pan Cytokeratin [Lu-5]
CM/PM/IP/ OAI162	Pan Cytokeratin Plus [AE1/AE3+ 8/18]
ACI/API/OIA/AVI3227	Pan Cytokeratin Plus [AE1/AE3+ 8/18]
VLTM162	Pan Cytokeratin Plus [AE1/AE3+ 8/18]
OPAI162	Pan Cytokeratin Plus [AE1/AE3+8/18]
AVI162	Pan Cytokeratin Plus [AE1/AE3+8/18]
API3035	Pan Lymphoma Cocktail (LCA+CD20+CD3+CD43)
PM362DS	Pan Melanoma + Ki-67
OPAI362	Pan Melanoma + Ki-67
PPM213DS	Pan Melanoma + S100 (Tyrosinase + MART-1 + S100)
CM/PM/OAI178	Pan Melanoma Cocktail -2 [Tyrosinase and Mart-1]
VLTM178	Pan Melanoma-2 Cocktail
S2801	Pascal Quality Strips
OPAI207	PAX5
VLTM207	PAX5
CM/PM/OAI207	PAX5
CP/PP379	PAX8
VLTR379	PAX8
OPAI379	PAX8
ACI/API/ OAI/AVI438	PAX8 (M)
VLTM438	PAX8 (M)
OPAI438	PAX8 (M)
PBS940	PBS Plus, 10X
VLTM3137	PD-1
ACI/API3137	PD-1
API3162	PD-1 (RM)



Biocare Medical LLC

ACI/API3224	PD-1 [CAL20]
ACI/API3171	PD-L1
OPAI3171	PD-L1
AVI3171	PD-L1
VLTR3171	PD-L1 (RB)
PXA969	PeroxAbolish
PX968	Peroxidazed 1
PX970	Peroxidazed Diluent
ACI/API3130	pHH3 (RM)
API300	PHOSPHO-EGFR
CP/PP404	Phospho-Histone H3
CRM/PRM350	Placental Alkaline Phosphatase (PLAP)
VLTM344	PMS2
CM/PM/IPI/AVI/ OAI344	PMS2
NC499	Polymer Negative Control Serum (Mouse & Rabbit)
ACI/API/AVI3252	PRAME [EPR20330]
ORI6024	Probe Enhancer (RNA)
ACI/API/OAI424	Progesterone Receptor (PR) [16]
API343	Progesterone Receptor (PR) [PgR636]
ACI/API302	Progesterone Receptor (PR) [SP2]
PRB957	Pronase Buffer
AVI3206	Prostate Cocktail (CK5 + CK14 + p63)
PM364	Prostate Cocktail-2X (CK5+CK14+p63)
ACI/API3163	Prostein
PR960	Protease XXIV
CM/PM278	PTEN (Tumor Suppressor)
OPAI278	PTEN (Tumor Suppressor)
CM309	PU.1
ORI6043	Rabbit AP Detection
OPRI6043	Rabbit AP Detection
ORI6008	Rabbit HRP Detection
OPRI6008	Rabbit HRP Detection
PD905	Renaissance Background Reducing Diluent
PD904	Renoir Red Diluent
RV1000X	Reveal Decloaker, 10X
RV1000	Reveal Decloaker, RTU
PD912	Revival Series Sampler



Biocare Medical LLC

RI0213	RISH™ AP Detection Kit
RI0208	RISH™ Dual Detection Kit
RPI0001	RISH™ Epstein-Barr Encoded RNA (EBER) Probe
RI0207	RISH™ HRP Detection Kit
RPI0004	RISH™ Kappa Light Chain DNA Probe
RPI0005	RISH™ Lambda Light Chain DNA Probe
RI0209	RISH™ Retrieval, 10X
RAEC810	Romulin AEC Chromogen Kit
ACI/API3208	ROR gamma T
ACI/API/VLTR3240	ROS1 [EPMGHR2]
HTEP	Rubens Eosin-Phloxine
VLTM128	S100 [15E2E2]
CP/PP/OAI021	S100 Protein (P)
CM/PM128	S100 Protein [15E2E2] (M)
ACI/API3237	S100 Protein [4C4.9] (M)
AVI3237	S100 Protein [4C4.9] (M)
OPAI3237	S100 Protein [4C4.9] (M)
VLTM3237	S100 Protein [4C4.9] (M)
CM/PM/IP/OAI089	S100 Protein Cocktail
VLTR3010	S100P
ACI/API3010	S100P
AVI384	SALL4
CM/PM/OAI384	SALL4
ACI/API3225	SATB2 [EP281]
SBT	Slide Brite
ACI/API3242	SMAD4 [EP618Y]
VLTM001	Smooth Muscle Actin (SMA)
CM/PM/IP/OAI001	Smooth Muscle Actin (SMA)
CM/PM420	Smooth Muscle Myosin Heavy Chain
CM/PM372	Smoothelin
VLTM372	Smoothelin
VLTM3099	SOX10 (M)
ACI/API/OAI/IP/AVI 3099	SOX10 (M)
VLTR3099	SOX10 (M)
OPAI3099	SOX10 (M)
VLTM3120	SOX-11
AVI3120	SOX11 (M)
ACI/API3120	SOX11 (M)
ACI/API3109	SOX2



Biocare Medical LLC

BRI/OR4039	SSC Wash Buffer
CM/PM406	β-Catenin
STUHRP700	Starr Trek Universal HRP Detection System
ACI/API3244	STAT6 [YE361]
CM/PM275	Surfactant Apoprotein-A
AVI371	Synaptophysin
OPAI371	Synaptophysin
VLTM371	Synaptophysin
CM/PM/IP/OAI371	Synaptophysin
NM-HEM	Tacha's Automated Hematoxylin
HTBLU	Tacha's Bluing Solution
HTBLUX	Tacha's Bluing Solution, 10X
TWB945	TBS Automation Wash Buffer, 20X
TWB946	TBS Automation Wash Buffer, 40X
CP/PP134	Terminal Deoxynucleotidyl Transferase (TdT {P})
VLTM022	Thyroglobulin Cocktail
CM/PM022	Thyroglobulin Cocktail
CM/PM412	Thyroid Stimulating Hormone (TSH)
AVI130	TIA-1
OPAI130	TIA-1
CM/PM130	TIA-1
VLTR3243	TIGIT (RM) [BLR047F]
ACI/API3254	TIGIT [BC41]
ACI/API3045	Topoisomerase II alpha
ACI/API135	Treponema pallidum (Spirochete)
PP416	TRIM29
TBS942	Tris Buffer Saline (TBS) Plus, 10X
CM/PM/IP087	TTF-1
ACI/API/OAI3126	TTF-1 [SPT24]
VLTM3126	TTF-1 [SPT24]
OPAI3126	TTF-1 [SPT24]
PM425DS	TTF-1 + CK5
PPM / IPI 394DS	TTF-1 + Napsin A
API3078DS	TTF-1 + Napsin A (RM)
VLTM3239	TTF-1 + Napsin A (VALENT)
API3141DS	TTF-1 + p40 (cRM)
ORI6036K	Tubing Cleaning Kit
OPRI6036K	Tubing Cleaning Kit



Biocare Medical LLC

CM/PM002	Tumor Associated Glycoprotein [B72.3]
CM/PM/OAI155	Tyrosinase
ORI6013	Universal Negative Control Serum
NC498	Universal Negative Control Serum
OPRI6013	Universal Negative Control Serum
API3001DS	Uro-2™ (CK20 + p53)
PM370TS	URO-3 Triple Stain (CD44+p53) with CK20
OPAI3051	Uroplakin II
ACI/API/OAI/ AVI3051	Uroplakin II
API3094	Uroplakin II + Uroplakin III
VLTM3094	Uroplakin II + Uroplakin III
ACI/API3023	Uroplakin III
VLT8012	Val Aqua Rinse, 20X
VLT8003	Val AR-Hi, 5X
VLT8004	Val AR-Lo pH, 5X
VLT8011	Val Background Block
VLT8028	Val Clean
VLT8023K	Val DAB Away Kit
VLT8010	Val DAB Buffer
VLT8009	Val DAB Chromogen
VLT8008K	Val DAB Chromogen Kit
VLT8001	Val DePar
VLT8021	Val Fast Red Buffer
VLT8020	Val Fast Red Chromogen
VLT8019	Val Fast Red Chromogen kit
VLT8014	Val Hematoxylin
VLT8006	Val Mouse Secondary HRP
VLT8002	Val Peroxidase Block
VLT8035	Val Plex 1
VLT8022	Val Plex 2
VLT8027	Val Prep
VLT8005K	Val Tect HRP Dection Kit
VLT8029	Val Universal HRP Linker
VLT8007	Val Universal HRP Polymer
VLT8015	Val Universal Negative Control
VLT8013	Val Wash Buffer, 20X
VLT8095K	Val Zyme Pronase Kit
VLT8081	Val Zyme Trypsin Buffer



Biocare Medical LLC

VLT8082	Val Zyme Trypsin Concentrate
VLT8080K	Val Zyme Trypsin Kit
VLT8101 with VLT8102/VLT8103	Valent Automated Slide Stainer
VLT8037	Val-Tect Mouse AP Polymer
VLT8047	Val-Tect Rabbit AP Polymer
VLT8031K	Val-Tect Universal HRP Detection Kit
PD902	Van Gogh Yellow Diluent
BRI4001	V-Blocker
ACI356	VEGF
AVI048	Vimentin
CM/PM/IP/OAI048	Vimentin
CRM/PRM312	Vimentin
VPD901	VP Monet Blue Diluent
FR805	Vulcan Fast Red Chromogen Kit 2
WR806	Warp Red Chromogen Kit
WR806MCK	Warp Red Chromogen Kit (Custom)
ORI6012	Wash Buffer
OPRI6012	Wash Buffer
ACI3238	WT1
API3238	WT1
CM/PM258	WT1 (Wilms' Tumor)
ACI/API3238	WT1 [rWT1/857]
ACI259	ZAP-70 (LR)
FRI7323	ZipFISH™ Hybridization Buffer

End of Product List